

## Redox control and measurement in low-temperature (<450 °C) hydrothermal experiments

JING FANG<sup>1,2</sup> AND I-MING CHOU<sup>1,\*</sup>

<sup>1</sup>CAS Key Laboratory of Experimental Study under Deep-sea Extreme Conditions, Institute of Deep-sea Science and Engineering, Chinese Academy of Sciences, Sanya, Hainan 572000, China

<sup>2</sup>University of Chinese Academy of Sciences, Beijing 100871, China

### ABSTRACT

Redox control in hydrothermal experiments is routinely achieved through double-capsule and Shaw membrane techniques. These techniques control oxygen fugacity ( $f_{\text{O}_2}$ ) by imposing a defined hydrogen fugacity ( $f_{\text{H}_2}$ ) on a studied sample enclosed, together with  $\text{H}_2\text{O}$ , in a hydrogen membrane capsule made of Pt or Ag-Pd alloys. However, due to the low permeability of these membranes to  $\text{H}_2$  at low temperatures ( $T$ ), these techniques do not work efficiently below 450 °C. Here, we tested fused silica as a new hydrogen membrane and successfully applied it to monitor and control the redox states of studied samples at  $T$  down to 200 °C in hydrothermal experiments. Our results showed that 3, 8, 16, 36, 96, and 216 h are sufficient for a fused silica capillary capsule (FSCC) to reach osmotic equilibrium with the externally imposed 1 bar of  $\text{H}_2$  at 350, 300, 250, 200, 150, and 100 °C, respectively, and  $\text{H}_2$  pressures inside a FSCC was very close to the externally imposed values after osmotic equilibrium. By using FSCC as a hydrogen fugacity sensor, equilibrium  $\text{H}_2$  pressures for Ni-NiO- $\text{H}_2\text{O}$  and Co-CoO- $\text{H}_2\text{O}$  redox buffer assemblages at 250–400 °C and 1000 bar total pressure were measured. The equilibrated  $f_{\text{O}_2}$  calculated are consistent with those derived from previous literature. Besides, FSCC can be used as a sample container, where  $f_{\text{H}_2}$  and  $f_{\text{O}_2}$  of enclosed samples can be continuously controlled. Furthermore, FSCC is an ideal container for sulfur-bearing samples, and its transparency allows spectroscopic analyses of the sample. Our work extended the low- $T$  limit of previously well-developed redox control techniques and may open up a new research avenue in low- $T$  hydrothermal experiments.

**Keywords:** Redox buffer and control techniques, hydrothermal experiments, redox reactions, hydrogen fugacity sensor, fused silica capillary